

## Appendix 9 – Programming and Setting the Function Generator

### A9-1 The Function Generator

A Tektronix AFG310 arbitrary function generator set-up as shown in Figure A9-1 is used to generate trigger signals for the Redlake ES1.0 camera by delaying the laser front panel sync out trigger signal.




Figure A9-1 Function generator as used to trigger the camera.

Among the tricky and useful things to know are:

- The manual trigger button which allows to trigger the camera manually (excellent for testing without the laser running).
- CH1 must be activated (LED shining) for the waveform to be generated.
  - Press the button beneath the CH1 LED to activate the waveform (bright LED).

## A9-2 Setting the Waveform and Frequency

This guide assumes that the laser Q-switch trigger is running at 10Hz and that the appropriate waveform (Chapter A9-3) has been programmed to User1. A summary of the necessary settings is shown in Table A9-1.

[ FUNC ] USR1	[ FREQUENCY (Hz) ] <sup>1</sup> 10.19		[ AMPL (V) ] 5.000
[ OFFSET (V) ] 0.000	[ MODE ] TRIG	[ MODUL ] Off	[ PHASE (DEG) ] 0

Table A9-1 Function Generator Settings

- 1) Turn on the frequency generator.
- 2) Set the waveform function to User 1:
  - Press [FUNC]; use the up and down arrows [  $\hat{\wedge}$  ] [  $\hat{\vee}$  ] to select "USR1";
  - Press [ENTER].
- 3) Set the frequency to 10.19 Hz by:
  - Press [FREQ];
  - Press [Hz/s/V];
  - Set to "10.19" by using the keypad and press [ENTER].
- 4) Set the amplitude to 5V (TTL high level):
  - Press [AMPL];
  - Set to "5.000" by using the keypad and press [ENTER]
- 5) Check if offset is zero, if not:
  - Press [OFFSET]; set to "0.000" as above.
- 6) Set to triggered mode
  - Press [MODE];
  - Set to "TRIG" by using the up and down buttons [  $\hat{\wedge}$  ] [  $\hat{\vee}$  ];
  - Press [ENTER].
- 7) Make sure that MODUL is OFF and that the PHASE is 0
  - If not correct in similar fashion as above.
- 8) Press CH1 (green light must appear) for operation!!!

<sup>1</sup> It might become necessary to change the frequency should any laser timing settings change. In this case the trigger signal delay has to be measured again (using an oscilloscope and a fast photodiode).

### A9-3 Programming the Waveform

The following guide shows how to program the waveform shown in Figure A9-2 (used to trigger the camera).




Figure A9-2 Waveform as programmed.

- 1) Set the waveform function to Edit:
  - Press [FUNC]; use the up and down arrows [  $\wedge$  ] [  $\vee$  ] to select "EDIT";
  - Press [ENTER].
- 2) Set the desired number of points:
  - Press [SHIFT], [EDIT] (same button as offset, blue letters);
  - Select 'NUM OF POINTS' by using the left and right buttons [  $\ll$  ] [  $\gg$  ];
  - Press [ENTER];
  - set value to "9990" by pressing [9], [9], [9], [0], and [ENTER]
- 3) Draw the waveform
  - {lines between points (1,2047) (9988,2047) (9988,4094) (9990,4094)}
  - Press [CANCEL] to move the cursor to 'NUM OF POINTS' again;
  - Use [  $\ll$  ] and [  $\gg$  ] to select 'LINE';
  - Press [ENTER];
  - Type [1], [ENTER];
  - Type [2], [0], [4], [7], [ENTER] → (1,2047)
  - Repeat for (9988, 2047)
    - (i.e. type [9], [9], [8], [8], [ENTER]; type [2], [0], [4], [7], [ENTER]);
  - and (9988,4094), and (9990,4094).

## Appendix 9 – Programming and Setting the Function Generator

### 4) Save the waveform to USR1:

- Press [CANCEL] to move the cursor to 'LINE' again;
- Use the left and right buttons [«] [»] to select 'SAVE TO';
- Press [ENTER], select 'USER1' by using the up and down buttons [ ^ ] [ v ]
- Press [ENTER] (this can take a short while).

### 5) Done!

- Press [CANCEL] (twice or until beeps occur) to get back to the main menu;
- And set up the parameters (frequency, .... → Chapter A9-2).